



Technical Specifications

MOBOTIX c71 Kepler NurseAssist *powered by Kepler Vision Technologies*

Detect falls and unauthorized bed leaving immediately in real time!

From now on, falls that are detected too late and unauthorized leaving of the bed can be avoided. Intelligent MOBOTIX sensor technology supports hospitals, nursing homes and retirement homes effectively and cost-efficiently in protecting patients and those in need of care - automatically around the clock.

- Monitor patient well-being with one optical sensor.
- Hemispherical optical sensor with 360° all-round vision without blind spots
- One sensor per patient/resident room is sufficient.
- Integrated AI-powered software detects incidents based on live images.
- Automatic alerts to staff in case of alarm.
- Direct patient dialog possible via integrated microphone/audio function.
- Simple and decent installation.
- GDPR compliant operation through event-controlled obscuring of live and recorded images.



Table of Contents

| | |
|--|-----------|
| Order Information | 2 |
| Hardware | 2 |
| General Software Features | 5 |
| App Specifications | 6 |
| Further Video Analysis | 8 |
| Video Management Software | 9 |
| Dimensions | 10 |

Order Information

| | |
|-------------|--|
| Name | MOBOTIX c71 Kepler NurseAssist powered by Kepler Vision Technologies |
| Order Code: | Mx-c71A-NurseAssist |

Hardware

| | |
|-----------------------------|---|
| Image sensor (DN) | 4k (1:1) effective max. resolution 2144x2144 |
| Light sensitivity | <ul style="list-style-type: none">▪ Color sensor (day): 0,1 lx @ 1/60 s; 0,005 lx @ 1 s▪ BW sensor (night): 0,02 lx @ 1/60 s; 0,001 lx @ 1 s |
| Exposure control | Manual and automatic mode 1 s to 1/16,000 s |
| Operating temperature range | -10 to 50 °C/-14 to 122 °F |
| Internal DVR Storage | Internal microSD card (SDHC/SDXC), 8 GB out of the box, max. 1 TB. |
| Microphone/Speaker | Integrated Microphone: <ul style="list-style-type: none">▪ Sensitivity: -35 dB +/-4 dB (0 dB = 1 V/pa, 1 kHz) Integrated Loudspeaker: |

| | |
|-----------------------------------|--|
| | <ul style="list-style-type: none"> ▪ 0.9 W at 8 Ohm |
| Integrated illumination | 850 nm IR, up to 15 m |
| Shock detector (tamper detection) | Yes |
| Max. power consumption | max. 25 W |
| Electrical surge protection | Available with MX-Overvoltage-Protection-Box (not part of the scope of delivery) |
| PoE standard | PoE Plus (802.3at-2009)/Class 4 (Class 3 can be enabled via HW switch. Note: IR illumination requires Class 4) |
| Interfaces | <ul style="list-style-type: none"> ▪ Ethernet 1000BaseT (RJ45 according to EIA/TIA-568B) ▪ USB-C/USB2.0; High-Speed ($V_{out} = 5.1\text{ V}$, $I_{out} = 0.9\text{ A}$, $P_{out} = 4.5\text{ W}$) |
| Mounting Options | Ceiling mountable |
| Dimensions (Ø x height) | 164 x 47,5 mm |
| Weight | Gross weight: 1.040 g |
| Housing | Aluminum, PBT-30GF |
| Standard accessories | <ul style="list-style-type: none"> 2 terminal blocks 1 mounting cable with sealing, 0.5 m/~2 in 2 silicone sealings for patch cable 1 Ethernet patch cable RJ45, 50 cm/19.7 in 1 Allen wrench 2.5 mm 1 c71 Kepler NurseAssist lens wrench 4 dowels 8 mm 4 stainless steel pan head-screws 4.5x60 mm 4 stainless steel washers Ø5.3 mm 1 Torx wrench TX 20 4 sealing plugs, white 1 Important Safety Information 1 Camera Connection Flyer |
| Detailed technical documentation | www.mobotix.com > Support > Download Center > Marketing & Documentation |

Technical Specifications

| | |
|-----------------------|--|
| MTBF | 80,000 hours |
| Certificates | EN 50121-4, EN 55032, EN 55035, EN 61000-6-1, EN 61000-6-2, EN 61000-6-3, EN 61000-6-4, EN 62368-1, EN 63000, AS/NZS CISPR32, 47 CFR Part 15b, NRTL |
| Protocols | DHCP (client and server), DNS, ICMP, IGMP v3, IPv4, IPv6, HTTP, HTTPS, FTP, FTPS, NFS, NTP (client and server), RTP, RTCP, RTSP, SFTP, SIP (client and server), SMB/CIFS, SNMP, SMTP, SSL/TLS 1.3, TCP, UDP, VLAN, VPN, Zero-conf/mDNS |
| Manufacturer warranty | 5 years |

General Software Features

| Feature | Properties |
|------------------------------|--|
| Software features | <ul style="list-style-type: none"> ▪ H.264, H.265 Multistreaming ▪ Multicast stream via RTSP ▪ Digital pan, tilt, zoom/vPTZ (up to 8x zoom) ▪ Genetec protocol integration ▪ Programmable exposure zones ▪ Snapshot recording (pre/post-alarm images) ▪ Continuous recording ▪ Event recording ▪ Time-controlled flexible event logic ▪ Weekly schedules for recordings and actions ▪ Event video and image transfer via FTP and email ▪ Playback and QuadView via web browser ▪ Animated logos on the image ▪ Master/Slave functionality ▪ Privacy zone scheduling ▪ Remote alarm notification (network message) ▪ Programming interface (HTTP-API) ▪ MxMessageSystem |
| ONVIF compatibility | Profile G, S, T, (M with later firmware release) |
| Master/Slave functionality | Yes |
| Remote alarm notification | email, network message (HTTP/HTTPS), SNMP, MxMessageSystem, MQTT |
| DVR/image storage management | <ul style="list-style-type: none"> ▪ On internal microSD card ▪ On external USB and NAS devices ▪ Different streams for live image and recording ▪ MxPEG+ only ▪ MxFFS with buffered archive, pre- and post-alarm images, storage monitoring with error reporting |

Technical Specifications

| Feature | Properties |
|---------------------------|--|
| Camera and data security | User and group management, SSL connections, IP-based access control, IEEE 802.1X, intrusion detection, digital image signature |
| Digitally signed firmware | Yes (to prevent firmware file tampering) |

App Specifications

General App Information

Product Name

Order Code

(only available in a bundle with the MOBOTIX c71 Kepler NurseAssist camera)

Camera & Basic App Bundle:

Mx-APP-KEP-NUA-EXT

Extended App license:

Mx-APP-KEP-NUA-EXT

Supported MOBOTIX Cameras

MOBOTIX c71 Kepler NurseAssist (basic version pre-installed)

Minimum Camera Firmware

v7.3.5.x

MxManagementCenter compatibility

- min. MxMC v2.5.3
- Configuration: Advanced Config license required
- Event Search: Smart Data Interface license included

General Features

| | |
|-----------------------------------|--|
| App Features Basic Version | <ul style="list-style-type: none"> ■ monitoring of the well-being of patients in private rooms, corridors and common areas of Healthcare facilities. ■ minimization of reaction times due to decentralized state-of-the-art fall detection. ■ prevention of accidents through person down, person in bed and person out of bed event. ■ automatized detection and localization of the bed in the patient's room. ■ GDPR compliant operation through event-controlled obscuring of live and recorded images. ■ Dynamic masking of a person in bed as additional privacy protection measure. |
| App Features Extended Version | <p>Plug-In license required to extend the NurseAssist status and alarm notifications with the following functions:</p> <ul style="list-style-type: none"> ■ Person sits on the edge of the bed ■ Person has left the room ■ Person has entered the room ■ Person has entered the bathroom / is in the bathroom ■ Bed is not recognized / is not present <p>NOTE! The Plug-In license requires MOBOTIX c71 NurseAssist Smart Sensor or licensed Kepler NurseAssist Certified App.</p> |
| Maximum number of detection areas | 1 |
| Recognized objects | Persons |
| Supported image sensor types | Day/Night MOBOTIX c71 Kepler NurseAssist default |
| Meta Data / Statistic format | JSON |
| Supported Reporting Interfaces | JSON / XML via HTTP(S), MQTT, MxMessageSystem |
| MxMessageSystem supported | Yes |

Technical Specifications

ported

| | |
|----------------|-----------------------------|
| MOBOTIX Events | Yes |
| ONVIF Events | Yes (Generic Message event) |

Scene Requirements

| | |
|--|---|
| Room size | max. 25m ² assuming an almost square room (e.g. 5x5 m) by radius of 4 m - long and narrow rooms will not work. |
| Required installation position (camera) | ceiling-mounted |
| Recommended installation height (camera) | 2,0m - 3,5m |
| Mounting distance from bed | 0,5m - 3,0m |
| Mounting distance from wall | > 70cm |

Technical App Specifications

| | |
|--------------------------------------|---|
| Synchronous / Asynchronous App | asynchronous |
| Simultaneous execution of other apps | No |
| Accuracy | min. 99% (considering scene requirement and no view blocking objects) |

Further Video Analysis

| Feature | Properties |
|------------------------|-----------------------------|
| Video motion detection | Yes |
| MxActivitySensor | Version 1.0, 2.1, 3.0 |
| MxAnalytics | Yes with a later FW release |
| MOBOTIX App support | Yes with a later FW release |



Video Management Software

| Feature | Properties |
|------------------------|---|
| MOBOTIX HUB | Yes www.mobotix.com > Support > Download Center > Software Downloads |
| MxManagementCenter | Yes (latest version recommended) www.mobotix.com > Support > Download Center > Software Downloads |
| MOBOTIX Cloud | Streaming & Events supported |
| MOBOTIX LIVE App | Yes (available in Google Play Store (Android) and Apple App Store (iOS)). |
| 3rd Party VMS Software | see ONVIF Profile S, T and G specification |

Dimensions

NOTE! Download the drilling template from the section or on the MOBOTIX website: www.mobotix.com > [Support](#) > [Download Center](#) > [Marketing & Documentation](#) > [Drilling Templates](#).

CAUTION! Always print or copy the drilling template at 100% of the original size!



MOBOTIX

BeyondHumanVision

EN_03/25

MOBOTIX AG • Kaiserstrasse • D-67722 Langmeil • Tel.: +49 6302 9816-103 • sales@mobotix.com • www.mobotix.com

MOBOTIX is a trademark of MOBOTIX AG registered in the European Union, the U.S.A., and in other countries. Subject to change without notice. MOBOTIX do not assume any liability for technical or editorial errors or omissions contained herein. All rights reserved. © MOBOTIX AG 2023